



4-year guarantee

PARI SINUS2

Technical datasheet

Precise, effective and gentle treatment of sinusitis.



Technical data:

Alternating current: 220 – 240 V, 50 Hz, **Pressure:** 1.5 bar, **Weight:** 1.7 kg, **Housing dimensions:** 18.5 x 13 x 15 cm (W x H x D), **Fill volume:** min. 2 ml, max. 8 ml, **Noise volume:** 55 db(A)

Second to none – pulsating aerosol effectively deposits the medicine or saline solution in the paranasal sinuses ^{1, 2}



Gamma camera images overlap with MRI images

The PARI SINUS2 includes:

- PARI LC SPRINT nebuliser



	NASAL SPRAY	PARI SINUS2
Nasal deposition in %	96.5*	21.5*
Deposition in the paranasal cavities in %	0.45*	12.5*

◀ Around 28 times more effective than the nasal spray

* Deposition (in %) of the applied active agent.

1. Möller W. et al (2009): Ventilation and aerosolized drug delivery to the paranasal sinuses using pulsating airflow – a preliminary study. Rhinology, 47, 405-412, 2009.
2. Schuschnig U. et al (2006): Comparison of delivery efficiency in a nasal case model of fluticasone propionate suspensions and a novel solution aerosolized via the PARI Vibrent™.

PARI SINUS2

Item numbers for your region

Product	Sales Region	IFU languages	REF (Item No.)
 PARI SINUS2 Nebuliser System	Western Europe	de, fr, it, pt, es, nl, en, da, fi, no, sv	128G1001
	Eastern Europe	cs, et, el, pl, ro, sl	128G1006
	UK & Ireland	en	128G1004
	Australia & New Zealand	en	128G1015
	Russia	ru	128G1018
 PARI SINUS2 Year Pack	Western Europe	de, fr, it, pt, es, nl, en, da, fi, no, sv	023G3811
	Eastern Europe	cs, et, el, pl, ro, sl	023G3816
	Russia	ru	023G3818

Useful accessories for PARI SINUS2

Product	REF (Item No.)
PPARI LC SPRINT SINUS Nebuliser with PARI LC SPRINT SINUS Application Set	023G2800
PARI LC SPRINT SINUS Application Set	023G2840
PARI LC SPRINT SINUS Nasal Joining Piece	041E0573
PARI SINUS Tube System	041G4570
PARI LC SINUS Nose Plug	041B0576
PARI LC SPRINT Nebuliser	023G1001
PARI Bag	085E1100



PARI GmbH
 Moosstraße 3, 82319 Starnberg, GERMANY
 E-Mail: export@pari.de
 Phone: +49 (0) 8151 279 220
www.pari.com